

Quality of closed chest compression on a manikin in ambulance vehicles and flying helicopters with a real time automated feedback

Havel C, Schreiber W, Trimmel H, Malzer R, Haugk M, Richling N, Riedmüller E, Sterz F, Herkner H.

Resuscitation

2010; 81(1):59-64

ARTICLE IDENTIFIERS

DOI: 10.1016/j.resuscitation.2009.10.007

PMID: 19926386

PMCID: not available

JOURNAL IDENTIFIERS

LCCN: not available

pISSN: 0300-9572

eISSN: 1873-1570

OCLC ID: not available

CONS ID: not available

US National Library of Medicine ID: not available

This article was identified from a query of the SafetyLit database.