

**Improved liquid chromatography-tandem mass spectrometric method for the determination of ethyl glucuronide concentrations in hair: Applications to forensic cases**

Imbert L, Gaulier JM, Dulaurent S, Morichon J, Bevalot F, Izac P, Lachâtre G.  
International journal of legal medicine  
2014; 128(1):53-58

**ARTICLE IDENTIFIERS**

DOI: 10.1007/s00414-013-0894-8  
PMID: 23824336  
PMCID: not available

**JOURNAL IDENTIFIERS**

LCCN: not available  
pISSN: 0937-9827  
eISSN: 1437-1596  
OCLC ID: 23022400  
CONS ID: sn 91033072  
US National Library of Medicine ID: 9101456

This article was identified from a query of the SafetyLit database.